ClarionMask

**CLASS I MEDICAL DEVICE
TYPE II SURGICAL MASK**

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ClarionMask

TECHNICAL SHEET

CLASS I MEDICAL DEVICE TYPE II SURGICAL MASK

COMPLIES WITH UNI EN 14683: 2014 AND UNI EN ISO 10993 STANDARDS (part 5 and part 10)

Medical Device Directory Number: 1965271

CLARIONMASK is a Class I and TYPE II Surgical Mask.

Like all TYPE II masks, it has a high filtering power in compliance with the reference standard EN 14683.

The traditional mask allows the deposit of viruses on the surface layers of the same, without eliminating them.

The superficial tissues of ClarionMask, on the other hand, are treated with LIVINGUARD Technology, the exclusive self-disinfecting formulation that mechanically inhibits the activity of bacteria and viruses and prevents their proliferation (eliminating them from the surface), guaranteed prolonged and effective protection from the risk of contamination.

ClarionMask is the mask that "inactivates" viruses,
including 99.9% of the SARS-CoV-2 which causes Coronavirus.

Technically, this technology applies a positive charge at the molecular level to the surfaces of tissues. When the microbes come into contact with these positive charges, the microbial cells, which have a negative charge, are eliminated, thus leading to the permanent destruction of the microorganisms.

PRODUCT DESCRIPTION

CLARIONMASK, 5-layer surgical mask with high filtering power by Clarion Living s.r.l. made in Italy, with exclusively Italian products and with a high sartorial quality, it consists of:

- an outer layer of colored cotton (100% poplin) fabric;
- 3 intermediate layers of polypropylene fabric according to an exclusive process called 3VR;
- an internal layer of colored cotton fabric (100% Jersey).

Each layer is structured differently for the specific function it is intended for: protective on the outside, filtering in the intermediate layers, protective and comfortable on the inside.

Inside the outer layer, in the upper central part, a nose shaping underwire is sewn, called "nose", in order to better adapt the closure of the mask in correspondence with the nose.

In the lower part, the adaptation to the face has been entrusted to an elastic webbing (cotton-based) inside which an elastic is sewn, the length of which allows it to pass behind the ear and return to be sewn in the upper side of the mask, in order to allow the perfect adherence of the mask to the face.

The elastics have a regulating ring for greater or lesser traction of the mask to the face, with two final rings that allow a further pull of the elastics behind the neck.

A fabric label is sewn inside the mask, which can be cut, showing the certifications of the same and its storage and washing methods.

On the outside and laterally on the left side, a cotton ribbon is sewn with the indications of the percentages of the products used for the realization of the CLARIONMASK and the logo of the Clarion Living manufacturer.

On the left side, the logo of the company that supplies the technology used for the antiviral treatment of fabrics (both external and internal), LIVINGUARD (Suisse), is hot stamped.

FABRIC TREATMENT

For the realization of CLARIONMASK, the highest quality Italian cotton fabrics were used. The same have been treated in Switzerland with Livinguard Technology.

Thanks to the use of these high quality fabrics, CLARIONMASK is extremely light, while at the same time ensuring a great filtration capacity for bacteria and viruses, together with a low resistance to respiratory flows.

All this makes the use of CLARIONMASK comfortable, while ensuring a high degree of safety, brought about by the permanent destruction of microorganisms.

CERTIFICATIONS – FABRIC TEST

The LIVINGUARD technology, applied to fabrics, renders ineffective up to 99.9% the viruses exhaled by humans, such as SARS-COV-2, known as CORONAVIRUS.

Freie Universität of Berlin (FUB), *RWTH Institute of Textile Technology* of Aachen and *University of Arizona* of Tucson with the tests carried out by them have unanimously certified that this disinfectant applied to fabrics is capable of definitively inactivating the Coronavirus that deposits on it.

REFERENCES

It is possible to consult on the LIVINGUARD or Clarion Living website the results of the official tests performed on:

- *Freie Universität* – Berlin
- Report test 30.07.2020 - Prof. Dr. Uwe Rösler and Dr. Anika Friese -
- *University of Arizona* – Tucson
- Report test 04.05.2020 - Luisa A. Ikner Ph.D. and Charles P. Gerba Ph.D. -

CERTIFICATIONS - CLARIONMASK

CLARIONMASK is a: **Class I Medical Device, Type II Surgical Mask**, according to the requirements of the UNI EN 14683: 2014 standard.


CLARIONMASK is registered with the Italian Ministry of Health, as a Medical Device (DM) at n ° 1965271, in the name of Clarion Living s.r.l.

USES

CLARIONMASK, is a reusable protective mask.

The outer fabric layer and the inner fabric layer are treated with the Swiss LIVINGUARD technology, whose technical characteristics are visible on the official website.

TECHNICAL FEATURES AND VALUES of ClarionMask

TECHNICAL FEATURES	VALUES																																
Sizes																																	
	Measurements in mm																																
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 10%;">SIZE</th> <th style="width: 10%;">A</th> <th style="width: 10%;">B</th> <th style="width: 10%;">C</th> <th style="width: 10%;">D</th> <th style="width: 10%;">E</th> <th style="width: 10%;">F</th> <th style="width: 10%;">G</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">L</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">M</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">S</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> </tbody> </table>	SIZE	A	B	C	D	E	F	G	L								M								S							
	SIZE	A	B	C	D	E	F	G																									
	L																																
M																																	
S																																	
Measurements in gr																																	
<table border="1" style="width: 100%; border-collapse: collapse;"> <tbody> <tr> <td style="width: 10%; text-align: center;">L</td> <td></td> </tr> <tr> <td style="text-align: center;">M</td> <td></td> </tr> <tr> <td style="text-align: center;">S</td> <td></td> </tr> </tbody> </table>	L		M		S																												
L																																	
M																																	
S																																	
Peso																																	

TECHNICAL FEATURES	VALUES	
1° colored outer layer	100% cotton POPLIN type fabric (Supima Cotton) - made in Italy -	
2°/3°/4° intermediate filter layer	3 superimposed layers of polypropylene fibrils - made in Italy -	
5° inner layer (in contact with the face)	100% cotton YERSEY type fabric - made in Italy -	
Edges	HIGHER	Tape sewn in 100% cotton - 20 mm
	INFERIOR	Tape sewn in cotton / elastane
Nose clip	Width 5 mm - iron rod Ø 0.5 mm - central plate 0.3 mm - covered with black polyethylene, conformable	
Nose clip length	mm 110	
Ear elastics	Empty tubular Ø 5 mm 83% Nylon and 17% Spandex - black hypoallergenic - elasticity 340% (± 5%) Length adjusters and final rings in Nylon and Silicon	
Tensile strength of the ties	Fm > N 35 = Kgf 3.57	
Sewing	Handcrafted in cotton	
Filtering efficiency (B.F.E.)	> 98%	Norm UNI EN 14683:2019+AC:2019 Lab. Techno Analisis – test report n. 2000579
Resistance to breath	< 18 Pa/cm2	Norm UNI EN 14683:2019+AC:2019 Lab. Techno Analisis – test report n. 2000579
Bioburden	< 11 UFC/g	Norm UNI EN 14683:2019+AC:2019 Lab. Techno Analisis – test report n. 2000579
Odours	Odourless totally	

TECHNICAL FEATURES	VALUES	
Other characteristics	<ul style="list-style-type: none"> • USE • FREQUENCY OF WASHING • MAXIMUM N° WASHING • EFFECTIVINESS • SHELF LIFE 	<p>> 8 hours a day 1 every 7 days 30 times 210 days 1 year from production date</p>
Substances which are not present	<ul style="list-style-type: none"> • Does not contain glass fibres • Does not contain latex • Does not contain PVC • Does not contain Phthalates • Does not contain Methylchloroisothiazolinone & Methylisothiazolinone • Does not contain dimethyl fumarate • Does not contain formaldehyde • Does not contain allergens • Does not contain heavy metals as: Antimony, Arsenic, Cadmium, Mercury and Lead • Nickel detected with method of artificial sweat: under limit fixed by Reach 	
Biocompatibility assessment	<ul style="list-style-type: none"> • Not irritant for skin - Norm UNI EN ISO 10993-10 – <i>Lab. Bio Basic Europe S.r.l.</i> • Does not show any contact dermatitis allergic potential - Norm UNI EN ISO 10993-10 – <i>Lab. Bio Basic Europe S.r.l.</i> 	
Cleanliness and Linting	<p>Cleanliness: 2,2 IPM Linting: 2,3 Log10</p>	
Sterility	<p>Not sterile</p>	

TECHNICAL FEATURES	VALUES
Primary packaging - Bag -	<ul style="list-style-type: none"> • Transparent bag • First protection of the mask • Bag made with 100% non-toxic and anti-suffocation polypropylene • Main technical data printed in two languages (ITA / EN) • Printing of the batch of origin, the production date and the expiry date • Adhesive closure
Secondary packaging - Case -	<ul style="list-style-type: none"> • Cardboard box with cream-colored background and Clarion Living green writing (for single sale) • Complete opening of the case so as not to damage the mask during collection • Man / woman model image print for better identification • Colored sticker for the identification of the color variant, to the main one • Printing of the batch of origin, the production date and the expiry date
Package leaflet	<ul style="list-style-type: none"> • Paper leaflet • Detailed information leaflet with technical characteristics and storage methods in two languages • Main technical data printed in six languages, of which two main (ITA / EN) and four secondary (SP / PT - FR / DE)
Characteristics of the production area	<p>CLARIONMASKs are produced in a laboratory that has obtained the certifications referred to in the ISO 9001: 2015 and ISO 13485: 2016 Standards</p>
Check on the final product	<p>During each production step</p>
Respect for the environment	<p>CLARIONMASK is produced with respect for the environment. The 58% of the mask is in natural fiber (100% cotton), while the remaining 42% is completely recyclable. In a year, each person uses at least 365 conventional masks. With the use of ClarionMask you need only 2.</p>
Standard of reference	<p>Class I - EU Regulation 2017/745 - Type II - UNI EN 14683: 2014 standard -</p>
Storage conditions	<p>Store in a healthy, cool and dry environment, away from heat sources, in odor-free environments with low humidity.</p>

Product validity	1 year from packaging
Washing method	<ul style="list-style-type: none">• by hand in running water• without any type of detergent• without wringing• do not dry in direct sunlight• do not use a dryer• do not iron

PACKING

Pieces in each	Case Carton Pallet	1 piece 180 pieces 4.320 pezzi	- 1 pcs x 24 cases 180 pcs x 24cartons
External dimensions	Case Carton Pallet	mm 120 x 18 x 120 mm 400 x 400 x 400 mm 800 x 1200 x 2000 (h)	
Grams of the carton	Case Carton	Cardboar Corrugated cardboard	GC1 - 325 g/mq KMFMT 24242 EB

TABLE OF THE CODES AND COLOURS IN WHICH CLARIONMASK IS AVAILABLE

REORDER CODE	COLOUR	PRIMARY PACKING (bag)	SECONDARY PACKING (case)
3CLG.M1188/3VR/clg32	BLUE	1 pcs	1 pcs
3CLG.M30363/3VR/clg32	BORDEAUX	1 pcs	1 pcs
3CLG.M500768/3VR/clg32	GREEN	1 pcs	1 pcs
3CLG.M4490/3VR/clg32	BLACK	1 pcs	1 pcs
3CLG.M700107/3VR/clg32	ORANGE	1 pcs	1 pcs

PRODUCTION PLACE

CLARIONMASK is made by Clarion Living and produced exclusively in the production laboratory which has obtained the following certifications for the DESIGN, PRODUCTION AND PACKAGING OF SURGICAL MASKS AND ...:

- QUALITY MANAGEMENT SYSTEM
Certificate N. :: IT52921A - ISO 9001: 2015;
- MEDICAL DEVICES - QUALITY MANAGEMENT SYSTEM
Certificate N. :: IT52921H - ISO 13485: 2016.

CUSTOMER CARE SERVICE

A customer assistance service is active and can be contacted for any information at the following address: info@clarion-living.com .

Comparison between **DPI FFP2** e **D.M CLASSE I TIPO II** masks (CLARIONMASK)

TEST	<u>FFP2</u> masks are devices for protection from solid particles	<u>D.M</u> medical devices are devices for protection against viruses and bacteria
-------------	---	--

	FFP2 MASK performance by legislation	CLARIONMASK tested performance
Filtering efficiency (B.F.E.)	≥ 95%	> 98%
Resistance to breath	≤ 60 PA/cm2	< 18 PA/cm2
Bioburden	-	Sì
Flow of exhalation air	95 L/MIN	> 180 L/MIN
Carbon dioxide evacuation	Sì	Sì

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Bearbeiter/in ur

30. Juli 2020

Report - Antiviral Activity of textile product HEFM47.CO.01

Tested by: Prof. Uwe Rösler and Dr. Anika Friese

Freie Universität Berlin, Institute for Animal Hygiene and Environmental
Health, Robert-von-Ostertag-Straße 7-13, 14163 Berlin

Tested textile product

- reference textile: Cotton (Co), untreated, 11.04.2020, lab number: C1
- antiviral textile: HEFM47.CO.01, 13.05.2020, LIVINGUARD, lab number: C2-neu

Method

- modified according to ISO 18184 (First edition 2014-09-01)
 - washing of all tested textiles with deionized water 10 times at 40°C, drying
 - cutting pieces of approximately 20x20mm and making up a mass of 0.4g with several pieces
 - sterilization at 121°C for 15 min, drying
 - before starting the test conditioning the textiles overnight in a humid environment (incubator 37°C)
 - controls: verification of cytotoxic effect and cell sensitivity to virus/inactivation of antiviral activity

- **test:** inoculation of 0,4g textile with 1ml (or 2ml) virus suspension (at least 10^7 TCID₅₀/ml)
 - washing out using 19 ml (or 18 ml) cell culture medium (MEM Eagle EBSS + 10% fetal calf serum + 1% non-essential amino acids, +1 % penicillin (10.000 U/ml) / streptomycin (10 mg/ml)) after specific time points
 - titration of washing solution in 96-well plate and titer calculation with method of Spearman and Karber in TCID₅₀/ml
- comparison of titer from untreated textile (reference) and treated textile for the specific time points

Results

Code ITU	C1	C2-neu	C1	C2-neu	C1	C2-neu	C1	C2-neu	C1	C2-neu	C1	C2-neu	C1	C2-neu
sequential nr.	1	1	2	2	3	3	4	4	5*	5*	6*	6*	6*	6*
textile after 30 min in log ₁₀ TCID ₅₀ /ml	7,6	6,35	7,6	6,6	7,35	6,35	7,35	6,475	7,35	7,1	7,475	6,975		
Antiviral Activity after 30min in log ₁₀		1,25		1		1		0,875		0,25		0,5		
Reduction after 30min in %		94,38		90,00		90,00		86,66		43,77		68,38		
textile after 1h in log ₁₀ TCID ₅₀ /ml	7,1	5,975	7,1	6,1	7,35	5,475	7,35	5,85	-	-	-	-		
Antiviral Activity after 1h in log ₁₀		1,125		1		1,875		1,5	-	-	-	-		
reduction after 1h in %		92,50		90,00		98,67		96,84	-	-	-	-		
textile after 2h in log ₁₀ TCID ₅₀ /ml	7,1	5,225	7,1	5,225	6,6	5,35	6,975	5,35	7,725	6,475	7,35	5,6		
Antiviral Activity after 2h in log ₁₀		1,875		1,875		1,25		1,625		1,25		1,75		
reduction after 2h in %		98,67		98,67		94,38		97,63		94,38		98,22		
textile after 6h in log ₁₀ TCID ₅₀ /ml	6,975	3,85	6,85	4,475	6,475	4,1	6,725	4,35	6,85	4,1	7,1	3,975		
Antiviral Activity after 6h in log ₁₀		3,125		2,375		2,375		2,375		2,75		3,125		
reduction after 6h in %		99,93		99,58		99,58		99,58		99,82		99,93		
textile after 12h in log ₁₀ TCID ₅₀ /ml	6,475	< 3,1	6,225	3,35	6,475	3,35	5,975	3,1	-	-	-	-		
Antiviral Activity after 12h in log ₁₀		>3,375		2,875		3,125		2,875	-	-	-	-		
reduction after 12h in %		>99,9		99,87		99,93		99,87	-	-	-	-		

*Inoculation with 2ml virus suspension instead of 1ml

Conclusion

The treated textile sample "C2-neu" can efficiently inactivate SARS-CoV-2. The results show a reduction of infectious SARS-CoV-2 within 6 hours from 2.3 to 3.1 log-ranges in comparison to the untreated textile. This means a reduction rate up to 99.9% after 6 hours. As a result of the inoculation of 0.4 g textile with 1 ml virus suspension the whole textile was moistened. This allows the textile sample having the largest disinfectant effect, we assume. Therefore, especially the handling with moist masks seems to be safer in comparison to the untreated textile.

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RESEARCH, DISCOVERY & INNOVATION

Water & Energy Sustainable Technology Center

Modified ISO 18184: Determination of Antiviral Activity of Textile Products to Evaluate Three Treated Fabrics against Human Coronavirus 229E

Company: Livinguard

Study Personnel:

Luisa A. Ikner, Ph.D.
Charles P. Gerba, Ph.D.

Test Date: 16 April 2020
Report Delivery Date: 04 May 2020

RESEARCH, DISCOVERY & INNOVATION

Water & Energy Sustainable Technology Center

Study Data

Table 1. Modified ISO 18184 Evaluation of Three Treated Fabrics against Human Coronavirus 229E at One Contact Time^{AM}

Test Virus	Contact Time	Sample ID	Virus Titer (TCID ₅₀ per Carrier)	Mean Virus Titer (TCID ₅₀ per Carrier)	Mean Log ₁₀ Virus Titer (TCID ₅₀ per Carrier)	Log ₁₀ Reduction	Percent Reduction
Human Coronavirus 229E (ATCC VR-740)	Time Zero	Control Fabric	2.25E+06	3.54E+06	6.55	N/A	N/A
			7.12E+06				
			1.26E+06				
	2 Hours	Control Fabric	7.12E+05	4.46E+05	5.65	0.90	87%
			4.00E+05				
			2.25E+05				
	2 Hours	FM-7 (Test Fabric)	2.25E+03	2.25E+03	3.35	2.30	99.5%
			2.25E+03				
			2.25E+03				
		FM-14 (Test Fabric)	4.00E+05	3.84E+05	5.58	0.08	13.8%
			4.00E+04				
			7.12E+05				
FM-26 (Test Fabric)		1.26E+04	7.34E+03	3.87	1.78	98.4%	
		7.12E+03					
		2.25E+03					

^{AM}TCID₅₀: Tissue Culture Infectivity Dose at the 50% Endpoint

^BLog₁₀ and Percent Reductions for Control Fabric at 2 hours calculated relative to Control Fabric immediately upon inoculation (Time Zero).

^CLog₁₀ and Percent Reductions for the three Test Fabrics at 2 hours calculated relative to Control Fabric mean viral titer of 2 hours.

RESEARCH, DISCOVERY & INNOVATION

Water & Energy Sustainable Technology Center

Table 2. ISO 18184 Evaluation of Three Treated Fabrics against Human Coronavirus 229E: Cytotoxicity Controls^{a,b,c}

Test Virus	Sample ID	Contact Time	Toxicity Titer (CCD ₅₀ per Replicate)	Mean Toxicity Titer (CCD ₅₀ per Replicate)	Mean Log ₁₀ Toxicity Titer (CCD ₅₀ per Carrier)
Human Coronavirus 229E (ATCC VR-740)	Control Fabric	2 Hours	1.26E+02	1.26E+02	2.10
			1.26E+02		
			1.26E+02		
	FM-7 (Test Fabric)	2 Hours	1.26E+02	1.26E+02	2.10
			1.26E+02		
			1.26E+02		
	FM-14 (Test Fabric)	2 Hours	1.26E+03	1.26E+03	3.10
			1.26E+03		
			1.26E+03		
	FM-26 (Test Fabric)	2 Hours	1.26E+02	1.26E+02	2.10
			1.26E+02		
			1.26E+02		

^aCCD₅₀: Cell Cytotoxicity Dose at the 50% Endpoint.

^bControl Fabric, FM-7 and FM-26 Test Fabrics: cytotoxicity observed for MRC-5 cells in the 10³ dilution.

^cFM-14 Test Fabric: cytotoxicity observed for MRC-5 cells in the 10⁵ and 10⁷ dilutions.

RESEARCH, DISCOVERY & INNOVATION

Water & Energy Sustainable Technology Center

Table 3. ISO 18184 Evaluation of Three Treated Fabrics against Human Coronavirus 229E: Neutralization Validation (NV) Controls^{a,b}

Test Virus	Sample ID	Contact Time	NV Ctrl Titer (TCID ₅₀ per Carrier)	Mean NV Ctrl Titer (TCID ₅₀ per Carrier)	Mean Log ₁₀ Virus Titer (TCID ₅₀ per Carrier)	Neutralization Validated?
Human Coronavirus 229E (ATCC VR-740)	Control Fabric	2 Hours	2.25E+03	2.72E+03	3.44	N.A.
			2.25E+03			
			4.00E+03			
	FM-7 (Test Fabric)	2 Hours	4.00E+03	3.30E+03	3.52	Yes
			2.25E+03			
			4.00E+03			
	FM-14 (Test Fabric)	2 Hours	4.00E+03	2.72E+03	3.44	Yes
			2.25E+03			
			2.25E+03			
	FM-26 (Test Fabric)	2 Hours	2.25E+03	2.25E+03	3.35	Yes
			2.25E+03			
			2.25E+03			

^aTCID₅₀: Tissue Culture Infectivity Dose at the 50% Endpoint.

^bNeutralization considered valid when mean viral titer on neutralized test fabrics differs by $\leq 0.5 \log_{10}$ relative to control fabric.

RESEARCH, DISCOVERY & INNOVATION

Water & Energy Sustainable Technology Center

Study Methods

Preparation of Control Fabric

1. An untreated, washed Control Fabric was provided by the Study Sponsor.
2. Pieces of the Control Fabric were cut to dimensions measuring 20 mm by 20 mm.
3. The Control Fabric cuttings were loaded into a glass beaker, covered with foil, and autoclaved for 20 minutes at 121 °C (103 kPa).

Preparation of Test Fabric

1. The treated Test Fabrics (FM-7, FM-14, and FM-26) were provided by the Study Sponsor.
2. Pieces of the Test Fabric were cut to dimensions measuring 20 mm by 20 mm.
3. The Test Fabric pieces were loaded into a separate glass beaker, covered with foil, and autoclaved for 20 minutes at 121 °C (103 kPa).

Test Procedure

1. On the day of testing, five to six 20 mm x 20 mm pieces each of the Control and Test Fabrics were aseptically transferred to and stacked within sterile Petri dishes. The final mass per container for the six swatches (Control or Test) was 0.40 ± 0.05 g. Nine Petri containers were prepared for the Control Fabrics, and six were prepared for each of the Test Fabric.
2. Six Control Fabric stacks and three Test Fabric stacks (per formulation) were each inoculated drop-by-drop with 0.2 mL of Human Coronavirus 229E viral stock (no soil load). A sterile pipette tip was used to press each stack and ensure that the inoculum was evenly spread through each piece of fabric.
3. Three Control Fabric and three Test Fabric stacks were parafilm and incubated at 20 °C in a humidified chamber for the 2-hour study contact time.
4. Three Control Fabric stacks were immediately neutralized to assess the viral titer upon inoculation (i.e. Time Zero) by transfer into conical tubes containing 4 mL of Lethen Broth. The tubes were then vortexed five times for five seconds each to wash out the viruses from the fabric pieces.
5. At the close of the 2-hour contact time, the triplicate Control and Test Fabric stacks were neutralized by transfer into conical tubes containing 4 mL of Lethen Broth. The tubes were vortexed five times for five seconds each to wash out the viruses from the fabric pieces.
6. The additional triplicate Control and Test Fabric stacks (containing no viral inoculum) were incubated concurrently for 2-hours and harvested in Lethen Broth as previously described to assess cytotoxicity and to validate neutralization as described in ISO 18184.

RESEARCH, DISCOVERY & INNOVATION

Water & Energy Sustainable Technology Center

Cell Culture Infectivity Assay

1. Control and test suspensions were diluted ten-fold in 0% FBS MEM.
2. Each dilution was plated in quadruplicate (0.1 mL per well) onto MRC-5 host cell monolayers in 24-well trays prepared to a confluency range of 70% to 80%.
3. Following an adsorption period of 30 minutes to facilitate virus-host cell interaction, 1 mL of 2% FBS MEM was added to each well.
4. The 24-well trays were incubated at 35 °C in a humidified chamber with an atmosphere of 5% CO₂ for 10 days. Assay trays were monitored regularly for changes to host cell monolayers indicative of cytotoxicity or viral cytopathogenic effects (CPE).
5. Assay trays were formally scored on Day 10 of incubation. Log₁₀ and percent reductions of the mean viral inoculum on the Test Fabrics at 2 hours were calculated relative to the mean viral titer of the Control Fabric at 2 hours using the Spearman-Kärber TCID₅₀ Method.
6. Log₁₀ and percent reductions were also determined for the Control Fabric at 2 hours relative to the titer yielded from Control Fabrics at Time Zero to ascertain stability of the viral inoculum over the course of testing.

Study Conclusions

The FM-7 Test Fabric was the most efficacious against Human Coronavirus 229E, achieving a reduction of 99.5% (2.30 log₁₀) given two hours of exposure relative to the Control Fabric (Table 1). Test Fabric FM-26 reduced levels of the test virus by 98.4%, while reductions were minimal for Test Fabric FM-14 at 13.8%.

With regard to cytotoxicity effects on the MRC-5 cell line, toxicity was observed in the 10⁰ dilution of washes for the Control Fabric, FM-7, and FM-26 (Table 2). Toxic effects on host cell monolayers were greater by an order of magnitude for FM-14 (10⁰ and 10⁻¹ dilutions). Neutralization was validated for each of the Control and Test Fabrics beyond the level of cytotoxicity, with < 0.5 log₁₀ difference in viral titer observed between the Control and Test fabrics (Table 3).

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San Felice sul Panaro (MO), 09/10/2020

TECHNO ANALISYS S.r.l.

Results are referred only to samples received in the laboratory.

This test report shall not be reproduced partially without written consent of the laboratory manager.

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CLASS I MEDICAL DEVICE

TYPE II SURGICAL MASK

BASED UPON THE TESTS CARRIED OUT IN COMPLIANCE WITH
UNI EN 14683:2019+AC:2019

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*All rights are reserved, being it a scientific and technical document protected by Copyright.
No part of this document may be reproduced in any form without the prior written authorization of Bio Basic Europe S.r.l.
Based on our experience, we suggest to check every 3 years its compliance with the guidelines in force.*

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ABSTRACT

To assess the potential skin irritation of the tested product, the sample underwent a process of extraction both in polar and non-polar solvent vehicles and the extracts were applied on reconstructed human epidermis (RHE) and then cell viability was assessed through a MTT test. The reduction of the viability of tissues exposed to the sample in comparison to negative controls is used to predict the skin irritation potential. Each extract was tested on 3 different RHE tissues. The test sample is considered to be irritant to skin if the tissue viability after exposure is $\leq 50\%$.

*The MTT assay showed that the viability of the tissues treated with the tested product is equal to **100.7%** for the extract in polar vehicle and **111.8%** for the extract in non-polar vehicle, showing that **tested product is not irritant for the skin.***

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INTRODUCTION

In vitro models to study human skin are important tools for research and development in the pharmaceutical and cosmetic industries. Human skin is the best possible model for such *in vitro* studies however, there are a number of legal and ethical issues concerning the use of human tissues. Animal skin is an alternative, but the relevance of conclusions drawn from animal data for human skin has always been questionable and moreover the EU regulation 1223/2009 prohibits use of animals for gathering toxicological data for cosmetic ingredients. In recent years several artificial human skin models have been developed: these models have undergone various testing in order to evaluate the possibility of using them to replace animal testing. For this purpose, they must mimic the relevant properties of human skin as closely as possible.

The RHE model consists of normal human keratinocytes cultured on an inert polycarbonate filter at the air-liquid interface. It represents a highly differentiated and stratified epidermis model comprising the main basal, supra basal, spinous and granular layers and a functional stratum corneum. The RHE model presents a histological morphology, composition and aspects of biochemistry comparable to the *in vivo* human tissue, so it is useful in several toxicity and permeability tests. Its use for skin irritation testing involves topical application of the product to the surface of the epidermis, and the subsequent assessment of its effects on cell viability.

Chemical-induced skin irritation, manifested mainly by erythema and oedema, is the result of a cascade of events beginning with penetration of the chemicals through the stratum corneum where they may damage the underlying layers of keratinocytes and other skin cells. The damaged cells may either release inflammatory mediators or induce an inflammatory cascade which also acts on the cells in the dermis, particularly the stromal and endothelial cells of the blood vessels. It is the dilation and increased permeability of the endothelial cells that produce the observed erythema and oedema. Notably, the RHE-based test methods, in the absence of any vascularisation in the *in vitro* test system, measure the initiating events in the cascade, e.g. cell / tissue damage, using cell viability as readout.

Cell viability is measured by enzymatic conversion of the vital dye MTT [3-(4,5-Dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide, Thiazolyl blue tetrazolium bromide], into a purple formazan salt that is quantitatively measured after extraction from tissues. Irritant chemicals are identified by their ability to decrease cell viability below defined threshold levels.

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AIM

The aim of the test was to evaluate the barrier effect of the tested product on reconstructed human epidermis (RHE).

The test was carried out on the liquid of maceration obtained incubating the product in saline (NaCl 0.9%) and sesam oil for 72 hours at 37°C.

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Tissues

Reconstructed Human Epidermis consists of normal human keratinocytes cultured on an inert polycarbonate filter at the air-liquid interface, in a chemically defined medium.

Trattamenti eseguiti / Treatments performed

NC	Controllo Negativo (D-PBS) / <i>Negative control (D-PBS)</i>
PC1	Controllo Positivo (SDS 1%) in solvente polare / <i>Positive Control (SDS 1%) in polar vehicle</i>
PC2	Controllo Positivo (SDS 1%) in solvente apolare <i>Positive Control (SDS 1%) in non-polar vehicle</i>
TS1	Estratto del campione in solvente polare <i>Extract of the sample in polar vehicle</i>
TS2	Estratto del campione in solvente apolare <i>Extract of the sample in non-polar vehicle</i>
V1	Solvente polare <i>Polar vehicle</i>
V1	Solvente apolare <i>Non-polar vehicle</i>

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Determination of cell viability by MTT assay

Cell viability was assessed by MTT assay. The chemical basis of the assay is the reduction of the MTT, a slightly yellow substance, which forms a purple formazan upon reduction. The process primarily takes place in the cytoplasm, and to a lesser extent in the mitochondria and cell membrane, and it is highly dependent on the concentration of intracellular NADH and NADPH. As a consequence of these metabolic processes, dark purple needle-like formazan crystals appear, radiating from the cells in a few hours. Formazan crystals can be solubilized in isopropanol. The absorption of the solubilized crystals can be measured at 540 nm wavelength and it is proportional to the viable cell number in an exceptionally wide linear range.

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Experimental protocol

Preparation of the extracts

The extracts were prepared according to ISO 10993-12 guidelines. We used as polar vehicle a saline solution of 0.9% NaCl, while as non-polar vehicle we used sesam oil. The extraction was performed at 37°C for 72 hours with continuous shaking.

Pre-incubation

RHE tissues were left in growth medium for at least 2 hours (37 °C, 5% CO₂) before treatment

Treatment

Each RHE tissue was treated with 100 µl of the extracts or of the controls. The tissues were then incubated at standard conditions for 18 hours. Each experiment was conducted in triplicate.

Rinsing

RHE tissues were repeatedly rinsed in D-PBS in order to remove all traces of the substances used for the treatments

MTT assay

RHE tissues were treated with MTT (1 mg/ml) for 3 hours and then the formazan crystals were extracted in isopropanol for 2 hours. Optical density read at 540 nm.

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The relative viability of each tissue was calculated on the basis of measured optical density (OD) values according to the following formula

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
The irritation potential is predicted if the relative viability of the tissues is ≤ 50% of NC

Trattamento <i>Treatment</i>	OD	Vitalità dei tessuti (%) <i>Tissue Viability (%)</i>
NC	1.329 ± 0.097	100.0 ± 7.32
PC1	0.065 ± 0.002	4.9 ± 0.18
PC2	0.073 ± 0.002	5.5 ± 0.15
TS1	1.338 ± 0.077	100.7 ± 5.77
TS2	1.486 ± 0.023	111.8 ± 1.70
V1	1.279 ± 0.036	96.2 ± 2.72
V2	1.412 ± 0.029	106.2 ± 2.17

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The values are expressed as means ± standard deviation deriving from the treatment of at least 3 tissues. NC = Negative control (SDS 1%); PC1 = Positive Control (SDS 1%) in polar vehicle; PC2 = Positive Control (SDS 1%) in non polar vehicle; TS1 = Extract of the sample in polar vehicle; TS2 = Extract of the sample in non polar vehicle; V1 = Polar vehicle; V2 = Non polar vehicle.

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ClarionMask 2 CLG M 1188/1S/clg-32 – taglia S

Is not skin irritant

The relative viability of the tissues treated with the extract in polar vehicle is 100.7% of negative control
The relative viability of the tissues treated with the extract in non-polar vehicle is 111.8% of negative control

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Kandarova H, Willoughby JA, De Jong WH, Letasiova S, Milasova T, Bachelor MA, Breyfogle B, Henda Y, De la Fonteyne L, Coleman KP. Pre-validation of an in vitro skin irritation test for medical devices using the reconstructed human tissue model EpiDerm™. *Toxicol In Vitro*. 2018 Aug;50:407-417.

Kandarová H. Dissertation: Evaluation and validation of reconstructed human skin models as alternatives to animal tests in regulatory toxicology. Free University of Berlin, 2006

Kandarová H, Liebsch M, Spielmann H, Genshow E, Schmidt E, Traue D, Guest R, Whittingham A, Warren N, Gamer AO, Remmele M, Kaufmann T, Wittmer E, De Wever B, Rosdy M. Assessment of the human epidermis model SkinEthic RHE for in vitro skin corrosion testing of chemical according to new OECD TG431. *Toxicology In Vitro*. 2006; 20: 547-59

ISO (2009) ISO/TC 194 10993-5 - Biological evaluation of medical devices – Part 5: Tests for in vitro cytotoxicity. International Organization for Standardization, Geneva, Switzerland.

ISO (2010) ISO/TC 194 10993-10 - Biological evaluation of medical devices – Part 10: Tests for irritation and skin sensitization. International Organization for Standardization, Geneva, Switzerland.

ISO (2012) ISO/TC 194 10993-12 - Biological evaluation of medical devices – Part 12: Sample preparation and reference materials. International Organization for Standardization, Geneva, Switzerland.

Kandarová H., Hayden, P., Klausner, M., Kubilus, J., Kearney, P. and Sheasgreen, J. (2009): In Vitro Skin Irritation Test: Improving the Sensitivity of the EpiDerm Skin Irritation Test Protocol. *ATLA* 37 (6), pp. 671 – 669.

Kandarová H., Willoughby, J.A., de Jong, W.H., Bachelor, M., Letasiova, S., Milasová, T., Breyfogle, B., de la Fonteyne, L., Coleman, K.P. (2015). Development, Optimization, and Standardization of an In Vitro Skin Irritation Test for Medical Devices Using the Reconstructed Human Tissue Model EpiDerm. *The Toxicologist*, Supplement to Toxicological Sciences 144, Abstract 2017, 2015.

Mosmann T. (1983). Rapid colorimetric assay for cellular growth and survival: application to proliferation and cytotoxicity assays. *J. Immunological Methods* 65, 55-62.

OECD (2015): OECD guideline for the testing of chemicals No 439: In Vitro Skin Irritation: Reconstructed Human Epidermis Test Method. Organisation for Economic Cooperation and Development, Paris, France.

DEF
Bio Basic Europe
Via G. Galvani, 10 - 20135 Milano
Tel. +39 02 57491111 - Fax +39 02 57491112
www.def.it

Bio Basic Europe
Via G. Galvani, 10 - 20135 Milano
Tel. +39 02 57491111 - Fax +39 02 57491112
www.bio-basic.com

Bio Basic Europe
Via G. Galvani, 10 - 20135 Milano
Tel. +39 02 57491111 - Fax +39 02 57491112
www.bio-basic.com

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Linee Guida OECD 442E

*In vitro analysis of the pro-sensitising potential
of a medical device*

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ABSTRACT

The pro-sensitising test was performed on cultures of monocytes treated with non cytotoxic concentration of the tested product. As positive control we used 2,4-dinitrochlorobenzene (DNCB), a well known sensitizing substance, while the negative control is represented by untreated cells maintained in culture medium. We evaluated the variation in the expression of two membrane marker, CD86 and CD54: a significant increase in the expression of these molecules on cells surface is indicative of the activation of an immune response and therefore the potential sensitizing action of the test substance.

In the cells treated with the tested sample we didn't observe a significant change in the expression of CD86 and CD54 compared to not treated cells. The lack of co-expression of the two markers means that the tested product has no pro-sensitising activity.

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INTRODUCTION

Skin sensitisation is the toxicological endpoint associated with chemicals that have the intrinsic ability to cause skin allergy in humans, resulting in a disease called Allergic Contact Dermatitis (ACD), a cell-mediated delayed-type hypersensitivity immune response induced by low-molecular-weight chemical compounds called haptens. The assessment of the skin sensitisation potential represents an important component in the safety evaluation of substances, especially those used in consumer products - e.g. cosmetics and detergents and represents a standard information requirement of legislation on chemicals in the EU. Traditionally, skin sensitisation hazard assessment has involved the use of laboratory animals. Within the European Union there is a prohibition of testing of cosmetic ingredients for the purposes of safety evaluation and the marketing of cosmetic products tested on animals is prohibited. Because of this regulation, and of the increasing social concern about animal welfare, the cosmetic industry has been developing alternatives to animal experimentation to test the efficacy or toxicity of cosmetic ingredients and other products.

The key biological events underlying the skin sensitisation process are well established: the activation of dendritic cells (DC), typically assessed by expression of cell surface markers, chemokines and cytokines, is considered to be a key event. DC can recognise and internalise antigens such as haptenated proteins, transport them via the lymphatic system to the regional lymph nodes and present them via major histocompatibility complex (MHC) molecules to naïve T lymphocytes to induce differentiation and proliferation of specific memory T-cells. In the lymph nodes, DC undergo a process of maturation characterised by phenotypic and functional changes resulting in the loss of their capability to process antigens and in the acquisition of the functionality of antigen presentation. The maturation process involves the decrease of phagocytic activity, increased expression of MHC molecules on the cell surface, changes in cytokines and chemokines secretion and up-regulation of several co-stimulatory (e.g. CD80, CD86, CD40) and intercellular adhesion molecules (e.g. CD54).

In the pro-sensitising test here proposed, the modulation of the CD86 and CD54 membrane markers in THP-1 cells, a human monocytic leukemia cell line used as a surrogate model for DC, is measured by flow cytometry following 24 hours of exposure to test sample. Significant variations in co-expression of these two markers with respect to not treated cells means that the tested sample has a pro-sensitising potential.

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AIM

*The aim of the test is to evaluate if the tested medical device does not cause pro-sensitising effects.
The test has been performed according to OECD Test Guidelines 442E.*

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Cell line and culture conditions

The test was carried out on human monocytes THP-1, used as surrogates for dendritic cells. It is known that THP-1 cells show enhanced CD86 and/or CD54 expression when treated with sensitizers. THP-1 cells were cultured in RPMI-1640 medium supplemented with fetal bovine serum (10%), 2-mercaptoethanol (0.05 mM) and antibiotics (1%). Cells were incubated at standard culture conditions (37°C, 5% CO₂). Good cell culture practices were used.

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Extraction fluid 50% and subsequent 1:2 dilutions

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Samples preparation

The tested product was incubated overnight in complete culture; the cells were treated with the liquid of maceration (filtered through a 0.22 µm) obtained in this way and subsequent dilutions 1:2. As positive control we used 2,4-dinitrochlorobenzene (DNCB) at the concentration of 4 µg/ml.

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Preliminary determination of cell viability

The product underwent a preliminary cytotoxicity screening on the THP-1 cells in order to choose a non-cytotoxic concentration to use in the pro-sensitising assay. The cell viability was evaluated through a MTS test. The MTS tetrazolium compound [3-(4,5-dimethylthiazol-2-yl)-5-(3-carboxymethoxyphenyl)-2-(4-sulfophenyl)-2H-tetrazolium, inner salt] is bioreduced by cells into a colored formazan product that is soluble in tissue culture medium. This conversion is presumably accomplished by NADPH or NADH produced by dehydrogenase enzymes in metabolically active cells. The assay was performed after a 24 hours-treatment period by adding a small amount of the MTS solution directly to culture wells, incubating for 1-4 hours and then recording the absorbance (optical density, OD) at 450 nm with a 96-well plate reader. The quantity of formazan product as measured by absorbance at 450 nm is directly proportional to the number of living cells in culture.

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Evaluation of the pro-sensitising potential

A suitable number of THP-1 cells was pre-cultured for 48 hours, seeded into a 24 well flat-bottom plate and treated for 24 hours with test sample and with the positive control. Not treated cells maintained in culture medium are negative controls. After 24 hours of exposure, cells were transferred into sample tubes and collected by centrifugation, supernatants were discarded and the remaining cells were washed twice with 600 µl of FACS buffer (PBS containing 0.1% BSA). Cells were split in two aliquots, centrifugated and resuspended in blocking solution (FACS buffer containing 0.01% globulin) and incubated at 4°C for 15 min. After centrifugation, cells were stained with FITC-labelled anti-CD86 or anti-CD54 at 4°C for 30 min. The expression levels of CD86 and CD54 were analysed using flow cytometry.

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Interpretation of results and prediction model

The expression of CD86 and CD54 is analysed with flow cytometry. Based on the geometric mean fluorescence intensity (MFI), the relative fluorescence intensity (RFI) of CD86 and CD54 for positive control (CTRL) cells and chemical-treated cells are calculated according to the following equation:

$$RFI = \frac{MFI \text{ of chemical-treated cells} - MFI \text{ of chemical-treated IgG1 cells}}{MFI \text{ of CTRL cells} - MFI \text{ of treated IgG1 CTRL cells}} \cdot 100$$

For CD86/CD54 expression measurement, each test chemical is tested in at least two independent runs to derive a single prediction (POSITIVE or NEGATIVE). The prediction is considered POSITIVE if at least one of the following conditions is met in 2 of 2 or in at least 2 of 3 independent runs, otherwise the prediction is considered NEGATIVE:

- *The RFI of CD86 is equal to or greater than 150% at any tested concentration (with cell viability \geq 50%);*
- *The RFI of CD54 is equal to or greater than 200% at any tested concentration (with cell viability \geq 50%).*

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The table shows the values for positive control (DNCB) and the 8 concentrations of the tested product

	DNCB	Concentrazione del campione / Sample concentration (%)							
		50	25	12,5	6,25	3,125	1,5625	0,7813	0,3906
CD86	152,16	111,64	107,33	112,50	103,02	105,60	106,47	101,72	106,03
CD54	204,76	107,14	109,52	111,90	107,94	109,52	112,70	115,87	118,25

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The calculated RFI values don't exceed the limits of the prediction model. We didn't observe co-expression of CD86/CD54

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CONCLUSIONS

The sample called
ClarionMask 2 CLG M 1188/1S/clg-32 – taglia S
has proved not to have pro-sensitising potential

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BIBLIOGRAFIA / REFERENCES

- Ainscough JS, Frank Gerberick G, Dearman RJ, Kimber I. Danger, intracellular signaling, and the orchestration of dendritic cell function in skin sensitization. *J Immunotoxicol* 2013; 10: 223-34
- An S, Kim S, Huh Y, Lee TR, Kim HK, Park KL, Eun HC. Expression of surface markers on the human monocytic leukaemia cell line, THP-1, as indicators for the sensitizing potential of chemicals. *Contact Dermatitis* 2009; 60(4): 185-92
- Ashikaga T, Yoshida Y, Hirota M, Yoneyama K, Itagaki H, Sakaguchi H, Miyazawa M, Ito Y, Suzuki H, Toyoda H. Development of an in vitro skin sensitization test using human cell lines: the human Cell Line Activation Test (h-CLAT). I. Optimization of the h-CLAT protocol. *Toxicol In Vitro* 2006; 20(5): 767-73
- Berridge MV, Tan AS. Characterization of the cellular reduction of 3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazolium bromide (MTT): Subcellular localization, substrate dependence, and involvement of mitochondrial electron transport in MTT reduction. *Arch. Biochem. Biophys* 1983; 303, 474-82
- Cory AH, Owen TC, Bartrop JA, Cory JG. Use of an aqueous soluble tetrazolium/formazan assay for cell growth assays in culture. *Cancer Commun* 1991; 3, 207-12
- Kaplan DH, Igyártó BZ, Gaspari AA. Early immune events in the induction of allergic contact dermatitis. *Nat Rev Immunol* 2012; 13: 114-24
- Kimber I, Basketter DA, Gerberick GF, Ryan CA, Dearman RJ. Chemical allergy: translating biology into hazard characterization. *Toxicol Sci* 2011; 1: 238-68
- Quah BJ, O'Neill HC. Maturation of function in dendritic cells for tolerance and immunity. *J Cell Mol Med* 2005; 9: 643-54
- Riss TL, Moravec RA. Comparison of MTT, XTT, and a novel tetrazolium compound for MTS for in vitro proliferation and chemosensitivity assays. *Mol Biol Cell* 1992 (Suppl.) 3, 184a
- Sakaguchi H, Ashikaga T, Miyazawa M, Kosaka N, Ito Y, Yoneyama K, Sono S, Itagaki H, Toyoda H, Suzuki H. The relationship between CD86/CD54 expression and THP-1 cell viability in an in vitro skin sensitization test--human cell line activation test (h-CLAT). *Cell Biol Toxicol.* 2009; 25(2): 109-26
- Sakaguchi H, Ashikaga T, Miyazawa M, Yoshida Y, Ito Y, Yoneyama K, Hirota M, Itagaki H, Toyoda H, Suzuki H. Development of an in vitro skin sensitization test using human cell lines; human Cell Line Activation Test (h-CLAT). II. An inter-laboratory study of the h-CLAT. *Toxicol In Vitro* 2006; 20(5): 774-84
- Tsuchiya S, Yamabe M, Yamaguchi Y, Kobayashi Y, Konno T, Tada K. Establishment and characterization of a human acute monocytic leukemia cell line (THP-1). *Int J Cancer* 1980; 26: 171-176
- Van Helden SF, van Leeuwen FN, Figdor CG. Human and murine model cell lines for dendritic cell biology evaluated. *Immunol Lett* 2008; 117: 191-197
- Vicini R, Mazzini G, Angolinotta C, Pintacuda V, Pastoris O. Forward scatter as a sensitive parameter to perform routine analysis of potential sensitizing compounds. XXXIII Conferenza Nazionale di Citometria, Lucca 22-25/09/2015
- Vocanson M, Hennino A, Rozières A, Poyet G, Nicolas JF. Effector and regulatory mechanisms in allergic contact dermatitis. *Allergy* 2009; 64: 1699-714

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RIPRODUZIONE RISERVATA © Copyright ANSA